

Mōveo™ XP mobile exercise platform

What is Mōveo XP?

The Mōveo XP is an innovative therapy platform for the medically complex or critically involved patient that allows for the early introduction of functional weight-bearing exercise in acute and post-acute care environments.

Why is early mobility an issue?

The devastating complications resulting from prolonged bed rest can affect all major organ systems, severely impact the patient's quality of life, and increase healthcare expenditures. A return to functional mobility is paramount for rehabilitation.

Why is Mōveo XP effective?

Deconditioning can be easily prevented when at-risk patients are treated early. Because partial weight-bearing exercise is ideal for this patient population, the Mōveo XP can become an invaluable tool to assist in the rehabilitation progression. As a mobile exercise platform, it is easily transported and utilized in the confines of an ICU, general patient room, and long term care centers. With a safe transfer method from bed or wheelchair, and customized mechanical patient support during exercise, the Mōveo XP allows for a transition to a weight-bearing stance in a safer, more comfortable and efficient manner. With the adjustability of the Mōveo XP, patients progress at their own rate while therapists document objective, incremental progress.



How does Mōveo XP work?

The Mōveo XP platform incorporates the benefits of tilt table standing with active exercise. Patients can access the platform directly from a wheelchair or from their hospital bed via a lateral transfer. The Mōveo XP's uniquely designed platform can be adjusted to an angle at which the patient is able to perform a controlled, shallow squat with comfortable resistance. As the patient's strength increases, the incline of the platform can be adjusted to provide more resistance until standing can be initiated with minimal assistance.

When should I use Mōveo XP?

Consider using the Mōveo XP for any patient cleared by a medical professional for mobilization, but who is limited by the following conditions:

- Lower extremity weakness
- Decreased LE ROM
- Decreased weight-bearing tolerance
- Poor postural control
- Poor endurance
- Poor muscle control

To learn more about the Mōveo XP, or to locate a dealer in your area, please call 800.592.7329, or visit chattgroup.com

